

PCTWORLD INTELLECTUAL PROPERTY ORGANIZATION
International Bureau

INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁵ : C12Q 1/54, G01N 33/00 A61K 37/00	A1	(11) International Publication Number: WO 92/07088 (43) International Publication Date: 30 April 1992 (30.04.92)
(21) International Application Number: PCT/US91/07629 (22) International Filing Date: 18 October 1991 (18.10.91) (30) Priority data: 600,797 22 October 1990 (22.10.90) US (71) Applicant: ABBOTT LABORATORIES [US/US]; CHAD 0377/AP6D-2, One Abbott Park Road, Abbott Park, IL 60064-3500 (US). (72) Inventors: SALDIVAR, Louis, Jr. ; 10606 Sheridan Road, Kenosha, WI 53143 (US). ENGLAND, Barbara, J. ; 1808 E. Kenwood Boulevard, Milwaukee, WI 53211 (US). (74) Agents: GORMAN, Edward, H., Jr. et al.; Abbott Laborat- ories, CHAD 0377/AP6D-2, One Abbott Park Road, Ab- bott Park, IL 60064-3500 (US).	(81) Designated States: AT (European patent), AU, BE (Euro- pean patent), CA, CH (European patent), DE (Euro- pean patent), DK (European patent), ES (European pa- tent), FR (European patent), GB (European patent), GR (European patent), IT (European patent), JP, KR, LU (European patent), NL (European patent), SE (Euro- pean patent). Published <i>With international search report.</i>	
(54) Title: STABILIZED ANOXIC DIAGNOSTIC REAGENT SOLUTION (57) Abstract A stabilized anoxic diagnostic reagent solution comprising an oxygen labile reagent, glucose oxidase, glucose, a hydrogen donor, sterile membrane fragments derived from bacteria having membranes containing an oxygen transfer system, and a reagent binding agent. The stabilized reagent has at least six months storage stability at 2-8 °C and open vial stability of at least three weeks.		

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AT	Austria	ES	Spain	MG	Madagascar
AU	Australia	FI	Finland	ML	Mali
BB	Barbados	FR	France	MN	Mongolia
BE	Belgium	GA	Gabon	MR	Mauritania
BF	Burkina Faso	GB	United Kingdom	MW	Malawi
BG	Bulgaria	GN	Guinea	NL	Netherlands
BJ	Benin	GR	Greece	NO	Norway
BR	Brazil	HU	Hungary	PL	Poland
CA	Canada	IT	Italy	RO	Romania
CF	Central African Republic	JP	Japan	SD	Sudan
CG	Congo	KP	Democratic People's Republic of Korea	SE	Sweden
CH	Switzerland			SN	Senegal
CI	Côte d'Ivoire	KR	Republic of Korea	SU ⁺	Soviet Union
CM	Cameroon	LI	Liechtenstein	TD	Chad
CS	Czechoslovakia	LK	Sri Lanka	TG	Togo
DE*	Germany	LJ	Luxembourg	US	United States of America
DK	Denmark	MC	Monaco		

⁺ Any designation of "SU" has effect in the Russian Federation. It is not yet known whether any such designation has effect in other States of the former Soviet Union.

STABILIZED ANOXIC DIAGNOSTIC REAGENT SOLUTIONField of Invention

The present invention is directed to a stabilized anoxic diagnostic reagent solution containing an oxygen labile reagent which may be stored in a liquid form at temperatures as high as 45° C for at least four (4) days, and permit their direct use as reference standards in instrumental analysis without the necessity for addition of water to them. In particular, the invention relates to the use of oxygen reducing enzyme systems in oxygen labile reagent solutions, especially reagent solutions used for calibrators and controls.

Description of the Related Art

Many materials deteriorate in the presence of oxygen. Various methods have been used to avoid the deteriorative effects of oxygen.

Biologically active substances such as found in sera, enzymes, lipids, hormones, electrolytes, and biologically active substrates or metabolites, are used widely in the diagnosis of diseases. They are used as reference standards for instrumental automated colorimetric analysis since they contain all or most of the components of the unknown to be analyzed. Once the diagnosing physician is aware of the basic concentration of components (differences vs. normal mean ranges of concentration of each components), the diagnosis can be made more objectively. In their natural form, when separated from their normal biological environment, such biologically active substances are unstable and undergo undesirable changes under the influence of heat, enzyme action, hydrolysis, oxidation and other influences causing undesirable molecular transformation therein. In the past, several methods of preservation have been utilized for such labile reagents.

One such procedure involves freeze-drying of the reagent. The freeze-drying procedure essentially involves rapidly reducing the temperature of the reagent containing solution followed by de-watering to a very substantial, if not total, extent at reduced pressure. If the reagent was additionally subjected to oxygen reduction then further procedures would be necessary. Not only is this procedure expensive, but it also requires reconstituting (resupplying) the normal water prior to use of the stored reagent. Natural proteins tend to

SUBSTITUTE SHEET

denature upon losing substantial concentrations of water. Thus, natural proteins could undergo compositional changes and loss of biological activity due to substantial de-watering occurring in the freeze-drying. Additionally, the freeze dried reagent often presents a turbid condition after it is reconstituted with water. This can cause analytical error because many of the automated analytical procedures basically involve colorimetry. During the reconstitution procedure, volume errors are sometimes introduced when, for example, the serum is reconstituted to its original liquid form by the addition of the water to the freeze-dried concentrate.

Prior to the discovery of this invention, certain diagnostic reagents such as bilirubin have been stabilized by polyols such as ethylene glycol. However, bilirubin stabilized in this manner must be kept at -20° C and it does not survive 45° C stress very well. Open vial stability of these materials is about 7 days. On the other hand, freeze-drying is the most frequently used method of stabilizing bilirubin. Once reconstituted, however, these materials do not have good stability, and errors in solubilizing the desiccated bilirubin may occur. These reconstituted standards generally require about 30 minutes for reconstituting.

U.S. Patent 4,414,334 to Hitzman discloses the removal of ambient oxygen from aqueous liquids, including canned foods and beverage products, which are catalyzed by an enzymatic deoxygenation system comprising alcohol oxidase in the presence of alcohol, optionally with the addition of catalase.

U.S. Patent 4,476,224 to Adler teaches the use of an enzyme system to promote the growth of anaerobic bacteria. The enzyme system of Adler comprises sterile membrane fragments derived from bacteria having membranes containing an electron transfer system which reduces oxygen to water. Optionally a hydrogen donor may be added to facilitate the reduction of water.

A technical bulletin, TB890411 released by Oxyrase, Inc. entitled "Properties of the Oxyrase™ Enzyme System" discloses that the Oxyrase™ enzyme system (sterile membrane fragments) in the presence of a suitable hydrogen donor, reduces dissolved oxygen directly to water. The bulletin also teaches that the enzyme system can be incorporated into oxygen-sensitive analytes used in automated, diagnostic instruments.

SUBSTITUTE SHEET

U.S. Patent 4,775,626 to Armenta, et al. teaches a method for reducing the oxygen content of a medium in which anaerobic cells are present. The method comprises having an effective amount of an oxidase and substrate for the oxidase in fluid contact with an aqueous medium. The medium can also contain a hydrogen peroxide scavenger.

The conversion of bilirubin through oxidation to biliverdin and other products is the main disadvantage of most bilirubin solutions. A system to prevent this oxidative process has been developed and is the invention presented here. The present invention offers a diagnostic reagent solution that is stable at 45°C for at least four days, has at least six months of open vial stability, and the finished product is already in liquid form.

SUMMARY OF THE INVENTION

The present invention is based upon the discovery of an aqueous bilirubin solution having sufficient stability that it can be shipped and stored for at least as long as six months at a temperature of 2-8°C and for at least four days at 45°C, has at least three weeks of open vial stability, and the finished product is already in liquid form.

→ The stabilization of the reagent is achieved by using a combination of
→ oxygen reducing enzyme systems in a synthetic matrix comprising a
→ reagent binding agent such as albumin to prevent the oxidation of the
→ reagent. The stabilized anoxic diagnostic reagent solution comprises an
→ oxygen labile reagent, a reagent binding agent, a first oxygen reducing
→ enzyme system which comprises glucose oxidase and its substrate
→ glucose, and a second oxygen reducing enzyme system which comprises
→ a hydrogen donor and sterile membrane fragments derived from
→ bacteria having membranes containing an oxygen transfer system
→ which reduces oxygen to water. The solution may also comprise
catalase, however, in amounts of 150 u/ml or greater open vial
stability decreases from that obtained with 50 u/ml. Additional
components, such as antioxidant used as a free radical scavenger,
alcohol, buffer and surfactant may also comprise the matrix. D-lactate
is the preferred hydrogen donor.

Over all, in the preferred embodiment bilirubin is stabilized by removing molecular oxygen, providing free, radical scavengers providing an anti-oxidant, providing a chelating agent, binding bilirubin to albumin, and providing for light protection in the final container. It is likely that several of these approaches work

SUBSTITUTE SHEET

synergistically with the enzyme systems to bring about the stability obtained. It is known that oxygen increases photodegradation, thus by removing oxygen we also reduce photodegradation. Reducing agents and chelators also behave synergistically with anti-oxidants. The advantages obtained are greater heat stability, the product is a liquid, the product has at least three weeks open vial stability, and the product does not contain human materials, thereby making this less hazardous than many of the commercially available standards which do contain human serum.

BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 is a graph comparing the activity of two biliruben A calibrators during 45°C heat stress testing.

Figure 2 is a graph comparing the activity of two biliruben B calibrators during 45°C heat stress testing.

Figure 3 is a graph comparing the activity of two biliruben C calibrators during 45°C heat stress testing.

THE DETAILED DESCRIPTION OF THE INVENTION

As used herein, and in the appended claims, the following terms, unless otherwise indicated, have the meanings set forth herein: "percent" and "parts" refer to percent and parts by weight; g small means gram or grams; mg means milligram or milligrams; l means liter or liters; ml means milliliter or milliliters; dl means deciliter or deciliters; ul means microliter or microliters; and M means molar and equals the number of moles of a solute in one liter of a solution; mM means millimolar. All temperatures are in degree Centigrade, °C.

The present invention involves the stabilization of a diagnostic reagent solution that comprises an oxygen labile reagent, or analyte, the addition of a reagent binding agent, and at least two oxygen reducing enzyme systems. It has been discovered that by adding an oxygen labile reagent to a matrix containing a reagent binding agent and the further addition of two oxygen reducing enzyme systems to that matrix a significant increase in the high temperature stability and open vial stability of the resulting reagent solution results.

In a preferred embodiment, a stabilized bilirubin solution that has been developed for the purpose of instrumentation calibration and

SUBSTITUTE SHEET

control is placed in a synthetic matrix. Enzyme systems are used to remove molecular oxygen. An antioxidant is used as a free radical scavenger and alcohol is used to solubilize the free radical scavenger wherein the alcohol may also be a known free radical scavenger. Albumin is used as the reagent binding agent to bind bilirubin which offers some additional stability.

The enzymatic reactions, of our invention, used for oxygen removal are as follows:

- 1) Glucose + Glucose oxidase + Oxygen - - >
Gluconic acid + Hydrogen peroxide
- 2) Hydrogen peroxide + catalase --- > oxygen +
water
- 3) Oxygen generated in step 2) recycles into step 1)
- 4) Membrane fragments + Oxygen + Hydrogen
doner - - > pyruvate + water

Inhibition is by the binding action of a reagent binding agent such as albumin, gluconic acid and possibly by a buffer system. The glucose, glucose oxidase reaction removes oxygen and prevents the formation of oxygen (O₂). Catalase may be used to remove hydrogen peroxide. The sterile membrane fragments, commercially available is Oxyrase, is a secondary oxygen removal system useful in the pre-reduction step of the following process.

By preventing reactive oxygen intermediates our invention should afford stability to many substances. Some substances which may have improved stability may be enzymes, reducing substances, pharmaceuticals, antibodies, antigens and substrates especially then linked to enzymes or other oxygen labile components.

EXAMPLES

The following examples, which described various modes including the best mode presently contemplated by the inventors, are presented solely for the purpose of disclosing and explaining the invention, and is not intended to be limiting.

EXAMPLE 1

The following is an example of the formulation process for producing a stabilized anoxic solution containing bilirubin as the oxygen labile reagent.

A matrix is prepared to contain the following in grams per liter.

a.	Ethanol	9.5 (to solubilize BHT)
b.	Polyoxyethylene 23 Lauryl Ether (BRIJ 35)	64.36
c.	Bistrispropane	30.037
d.	Glucose	5.3
e.	Butylated Hydroxy Toluene (BHT)	0.213
f.	FD-lactate	3.57
g.	Bovine Serum Albumin (BSA)	53.15
h.	Gentamicin	1.0
i.	HCL	As necessary to adjust pH
j.	NaOH	4 (To dissolve bilirubin powder)

The pH is adjusted to 7.85.

Oxyrase™ enzyme system (Oxyrase, Inc., Ashland, OH) is spiked into the matrix so as to obtain 0.3 units/ml of Oxyrase. (A unit of activity reduces 1.0% of the initially dissolved oxygen in each second per milliliter under standard conditions. The standard conditions are 37°C, pH of 8.4 in a 20 mM Phosphate buffer, 10 mM Sodium locate as the hydrogen donor, and an air saturated solution.) Oxyrase™ enzyme system is described in the literature and patents (4,476,224) as sterile membrane fragments derived from bacteria having membranes containing an electron transfer system (oxygen transfer) which reduces oxygen to water. Any such membrane fragments should operate similarly and could be substituted for Oxyrase. After the addition of the enzyme system room air must be excluded from the solution thereafter. This can be accomplished by maintaining nitrogen in the headspace.

Catalase from about 0 to about 150 units/ml can be added to remove any hydrogen peroxide which may be generated by the glucose oxidase reaction. However it has been observed that catalase in amounts more than 50 units/ml reduced open vial stability may result. After the oxygen has been removed by the Oxyrase™ enzymatic reaction or for about 60 minutes, glucose oxidase is spiked into the matrix so as to obtain 13 units/ml of glucose oxidase. Bilirubin dissolved in 0.1 M sodium hydroxide is then spiked in until the proper amount for that particular solution is achieved. The bilirubin powder

SUBSTITUTE SHEET

(Phansthiel, Inc., Waukegan, Illinois) dissolved by sodium hydroxide must not be hazy, but clear and amber colored prior to its addition to the matrix solution. The bilirubin may be conjugated, unconjugated or both.

As the following examples and data will show the following ranges for specific components added to a bilirubin calibrator solution are preferable in accordance with this invention.

1. pH - after all glucose is converted to gluconic acid (above 7.0)
2. Glucose - initial glucose level 0.1% to 0.5%, excessive glucose may cause bilirubin precipitation
3. BHT - 0% to about 0.024%
4. Oxyrase - 0.3 to about 0.6 units/ml
5. Albumin - 4% to 6%
6. Catalase - 0 to about 150 units/ml, excess catalase may reduce open vial stability.
7. Alcohol - about 1%
8. A non-ionic surfactant - about 5% to about 6%
9. Glucose oxidase - 10 to about 27 units/ml
10. Bilirubin - 0 to about 30 mg/dl
11. Hydrogen donor - 5 to about 20 mM

other ingredients may also be present. Table 1, below, shows the basic elements necessary to produce a stabilized anoxic reagent solution encompassed by our invention.

TABLE 1

<u>Ingredient</u>	<u>Amount</u>
Glucose	1-10,000mg/dl
Glucose oxidase	0.1-2 units/l
Binding agent	1-10%
Hydrogen donor	1-100 mM
Membrane fragment enzyme system	0.1-1 unit/l

Example 2

Bilirubin calibrators, at levels of 2 mg/dl(A), 10 mg/dl(B) and 20 mg/dl(C), were prepared by the method of Example 1 a set of vials of each calibrator level was kept at 45°C. Once each day the vials were tested and bilirubin levels were determined using a Vision™ analyzer and test cards. (Abbott Laboratories, Abbott Park, Illinois) according to the manufacture's instructions.

The control calibrators were stored at 2-8°C.

SUBSTITUTE SHEET

Figure 1, 2 and 3 show the results of the above study as well as the results of a parallel test conducted on a commercially available bilirubin calibrator (Medical Analysis Systems _____) prepared to comparable levels. In figure 1, 2 and 3 MAS CAL is the commercial calibrator and GO & OXY CAL is the calibrator of this invention. The points plotted on the graphs are average values. The calibrators prepared with the two enzyme systems of this invention showed superior stability and activity over the commercial formulation.

Example 3

Vials of bilirubin calibrators at each of three levels (as in Example 2) were prepared as in Example 1. Table 2 shows the exact composition of the solutions and the results of the open vial study, discussed below, in days. The reagent solution was considered acceptable if it was within the following ranges when the activity was compared to the appropriate control: calibrator A 15%, calibrator B 5%, calibrator C 5%

Three vials of each of the three calibrator levels were stored at 20- 80°C. Each work day (Monday through Friday) the vials were opened, then closed and left on a bench top for eight (8) hours. A baseline value was determined for each calibrator level by testing three vials from each level for five days prior to the commencement of this experiment. Throughout the test period bilirubin levels were obtained at random points in time.

Bilirubin values were obtained by using a Vision™ analyzer (Abbott Laboratories, Abbott Park, Illinois) and test cards according to the manufacture's instructions.


TABLE 2

Glucose	0.1%	0.3%	.5%
Glucose Oxidase	12.9 u/ml	2.9 u/ml	20 u/ml
Catalase	0	50 u/ml	150u/ml
Oxyrase	0.6 u/ml	0.3 u/ml	0.6 u/ml
DL Lactate	30 mM	30 mM	30 mM
BHT	0.02%	0.02%	0.02%
ETOH	1%	1%	0.95%
BRIJ 35	6%	5%	6.4%
pH	7.74	7.7	7.8

SUBSTITUTE SHEET

Bistrispropane	0.1 M	0.1 M	0.1 M
Gentamicin	0.1%	0.1%	0.1%
BSA	5%	5%	5%
Bilirubin	25 mg/dl	A,B,C	2A,2B,2C
Open Vial (Days)	40	25	15

In our formulation the following reactions are prevented. Therefore it is probable that many products may have extended stability in the Table 1 formulation.



The Haber-Weiss is prevented by the removal of molecular oxygen which prevents super-oxide formation. An alternate inhibitor of the Haber-Weiss reaction is the degradation of H_2O_2 . The Fenton reaction is inhibited by the chelation of iron.

The inhibition of these two reactions is believed to be a major source of stability in the formulation. To scavenge any preexisting free radicals BHT or BHA may also be incorporated into some formulations.

Example 4

Many complex biological molecules such as hormones, lipids, steroids, catechols, etc. as well enzymes, reducing substances, pharmaceuticals, antibodies, antigens and substrates could have their shelf and open vial life extended by this invention. The particular substance would place limits on the substrate concentrations of Table 1 and any buffer system which may be necessary.

A stabilized anoxic solution containing an oxygen labile reagent could be prepared as in Example 1 which would comprise the substrate of Table 1. Additionally the solution could further comprise BHT or BHA from 0.001 to 0.04%, buffer such as bistrispropane with adequate strength to maintain a pH compatible with the system, a detergent and organic solvent as required for BHT micel formation. Low levels of glucose oxidase are generally adequate since the substrate turnover is very high, however high levels may lead greater stability. The absorbance of this yellow compound may be a limiting factor in spectral applications.

SUBSTITUTE SHEET

Some enzymes and organic and inorganic chemicals are given stability by substances and reducing in nature. Cysteine and ascorbic acid, for example, when used in this fashion are of limited value since they are themselves subject to oxidation. In a matrix such as described in the above examples these oxygen labile substances could be stabilized so that they in turn could function as stabilizes.

In the case of enzyme reagents strong chelators may inactivate the enzyme by chelation of metallic coenzyme. In such a case the gluconic acid generation may be limited by limiting the glucose substrate. Degassing, pretreating with the membrane fragment enzyme system and then addition of glucose oxidase may also limit the generation of gluconic acid.

It will be appreciated by one skilled in the art that the embodiments described and the alternative embodiments presented are intended as examples rather than as limitations. Thus, the description of the invention is not intended to limit the invention to the particular embodiments disclosed, but it is intended to encompass all equivalents and subject matter within the spirit and scope of the invention as described above and as set forth in the following claims.

We Claim:

1. A stabilized anoxic diagnostic reagent solution that comprises:
 - an oxygen labile reagent;
 - a reagent binding agent;
 - a first oxygen reducing enzyme system wherein said first enzyme system comprises glucose and glucose oxidase; and
 - a second oxygen reducing enzyme system which comprises a hydrogen donor and sterile membrane fragments derived from bacteria having membranes containing an oxygen transfer system which reduces oxygen to water.
2. The anoxic reagent solution of claim 1 which further comprises catalase.
3. The anoxic reagent solution of claim 1 wherein the oxygen labile reagent is selected from the group consisting of hormones, lipids, steroids, catechols, pharmaceuticals, antibodies, antigens, substrates and reducing substances.
4. The anoxic reagent solution of claim 1 wherein the oxygen labile reagent is bilirubin.
5. The anoxic reagent solution of claim 1 wherein the reagent binding agent is albumin.
6. The anoxic reagent solution of claim 5 wherein albumin is bovine serum albumin.
7. The anoxic reagent solution of claim 1 that further comprises a free radical scavenger.
8. The anoxic reagent solution of claim 7 that further comprises an antioxidant.
9. The anoxic reagent solution of claim 8 wherein the antioxidant and the free radical scavenger comprises BAHT and ethanol.
10. The anoxic reagent solution of claim 1 wherein the hydrogen donor is d-lactate.

11. The anoxic reagent solution of claim 1 having a pH from about 7.4 to about 8.1.
12. A stabilized anoxic diagnostic reagent solution having a pH range from about 7.4 to about 8.1 and consisting essentially of a first oxygen reducing enzyme system, a second oxygen reducing enzyme system and bilirubin in a matrix of bistrispropane, glucose, d-lactate, NaOH, alcohol, BHT and bovine serum albumin, wherein said first enzyme system comprises glucose oxidase and catalase and said second enzyme system comprises sterile membrane fragments derived from bacteria having membranes containing an oxygen transfer system which reduces oxygen to water.
13. The anoxic reagent solution of claim 12 wherein the alcohol is ethanol.
14. The anoxic reagent solution of claim 13 wherein the ethanol is about 1 percent.
15. The anoxic reagent solution of claim 12 wherein the BAHT is present from about .0016 to about 0.024 percent.
16. The anoxic reagent solution of claim 12 wherein the bovine serum albumin is present from about 1 to about 6 percent.
17. The anoxic reagent solution of claim 12 wherein the alcohol is ethanol and is present from about 0.9 to about 1.1 percent.




FIG. 1




FIG. 2

SUBSTITUTE SHEET

2/2




FIG. 3

SUBSTITUTE SHEET

INTERNATIONAL SEARCH REPORT

International Application No. **PCT/US91/07629**

I. CLASSIFICATION OF SUBJECT MATTER (if several classification symbols apply, indicate all) *		
According to International Patent Classification (IPC) or to both National Classification and IPC IPC(5): C12 Q 1/54; G01N 33/00; A61K 37/00 USCL: 435/14,801; 436/97; 514/21		
II. FIELDS SEARCHED		
Minimum Documentation Searched ?		
Classification System	Classification Symbols	
US CL	435/14, 801; 436/97; 514/21	
Documentation Searched other than Minimum Documentation to the Extent that such Documents are Included in the Fields Searched *		
Chem Abstracts		
III. DOCUMENTS CONSIDERED TO BE RELEVANT *		
Category *	Citation of Document, ** with indication, where appropriate, of the relevant passages **	Relevant to Claim No. **
A	US, A. 4.414.334 (Hitzman) 08 November 1983. see entire document.	1-17
Y	US. A. 4.144.306 (Figueras) 13 March 1979. see column 16.	1-17
Y	US. A. 4.476.224 (Adler) 09 October 1984. see columns 1-3.	1-17
Y	US. A. 4.775.626 (Armenta et al.) 04 October 1988. see entire document.	1-17
Y,P	US. A. 5.047.395 (W) 10 September 1991. see column 1.	1-17
<p>* Special categories of cited documents: **</p> <p>"A" document defining the general state of the art which is not considered to be of particular relevance</p> <p>"E" earlier document but published on or after the international filing date</p> <p>"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)</p> <p>"O" document referring to an oral disclosure, use, exhibition or other means</p> <p>"P" document published prior to the international filing date but later than the priority date claimed</p> <p>"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention</p> <p>"X" document of particular relevance: the claimed invention cannot be considered novel or cannot be considered to involve an inventive step</p> <p>"Y" document of particular relevance: the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.</p> <p>"A" document member of the same patent family</p>		
IV. CERTIFICATION		
Date of the Actual Completion of the International Search	Date of Mailing of this International Search Report	
26 November 1991	30 JAN 1992	
International Searching Authority	Signature of Authorized Officer	
ISA/US	R. Girard	